

P R E C A U T I O N

Preserving

the safety

& dignity

P R E V E N T I O N

of patients

& caregivers

alike

P R O T E C T I O N

SWALTERNATE™

Alternation Therapy –

Provides alternation and low air loss therapies.



COMFORT

Change firmness of mattress with the touch of a button.

CARE

Two therapy modes: Static and 1 in 3 Alternation.

CONVENIENCE

35" width, 600 lb. capacity, and step deck system options.
(Bariatric sizes available.)



SIZEWISE

800.814.9389 | sizewise.net

SW ALTERNATE™

CONTROL UNIT

Liters of Air Per Minute	1740
Alternating Cycle Time	Programmable from 3 to 20 minutes
Static Low Air Loss	Yes
Pulsating Low Air Loss	N/A
Continuous Lateral Rotation	N/A
Alarms	Power Failure & Low Pressure
Warranty	2 year
Pressure Sensor Technology	Yes
Patented CPR Release	Yes
Dimensions	18" W x 7 1/2" D x 11 1/2" H
Unit Weight	16 lbs.


MATTRESS

Weight Capacity	600 lbs.
Mattress Size	Standard 35"
Percussion Modules	Add-on Feature
Warranty	1 year
Top Cover Material	200 Denier Polyurethane Coated
Height Inflated	7"
Mattress Weight	37 lbs.
Top Cover Warranty	90 days

ELECTRICAL

Rated Voltage	110 volts
Rated Frequency	50 Hz
Rated Input Power	120 volts
Fuse Rating	125 volts 5 amps
Power Failure Alarm	Yes

SW ALTERNATE™ CONTROL PANEL



800.814.9389 | sizewise.net