

BIGTURN™

Bariatric Rotation Therapy –
Provides full body lateral rotation with low air loss therapy.



PRECAUTION

Preserving

the safety

& dignity

PREVENTION

of patients

& caregivers

alike

PROTECTION

COMFORT

Change firmness of mattress with the touch of a button.

CARE

Adjustable modes and cycle times, side wall bolster support with < 40,> 40 degree rotation. Available with Percussion/Vibration option.

CONVENIENCE

Sizes to fit your needs: 39" or 48" widths. Accommodates weight capacities up to 1,000 lbs., scale available on the Bari-Rehab Platform™. (Non-bariatric sizes available.)




SIZEWISE

800.814.9389 | sizewise.net


SW BIG TURN II™

The Big Turn and the Big Turn II therapeutic mattress replacement systems are the only products on the market that provide full-body lateral rotation of 40 degrees or more for the bariatric patient weighting up to 1,000 pounds. The Big Turn and Big Turn II also provide the standard pressure redistribution therapy of low air loss to prevent and assist in the treatment of skin issues such as pressure ulcers. The newly enhanced Big Turn II provides the added pressure redistribution component of pulsation and an optional respiratory therapy module that allows mechanical percussion and vibration for the immobile patient weight up to 1,000 pounds.



PATIENT TURNING RIGHT ILLUSTRATION

PERCUSSION/VIBRATION



The purpose of chest physiotherapy (CPT) is to help patients breathe more freely and to facilitate oxygen and carbon dioxide exchange in the lungs. Percussion and vibration therapy is used to help break up thick secretions in the lungs by repeatedly stimulating the chest wall. Our percussion and vibration therapy can be programmed to focus on one, two or all three zones of the back, targeting the movement of fluid in the patient's lungs. The programmable therapy cycle can be set for five, ten, or twenty minutes with options for percussion, vibration, or both. This therapy is available on the Big Turn II or SW Rotate in 35", 39", or 48" widths.

MADE IN THE USA



SIZEWISE

800.814.9389 | sizewise.net